DISPOSITIVO MÉDICO

# PLEJUV

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# PLEJUV

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

[Reference]  
5) Chinese Chemical Letters 32(2021)577-582

inflammatory response shortly after the injection, while the inflammation was decreased or disappeared on Day 7 or Day 10 (Figs. 56 and 57 in Supporting information). Nevertheless, due to the rapid degradation rate of PDLLA MS, a large number of PDLLA MS were degraded and absorbed in the skin tissue on Day 30, which could not guarantee the long-term effect. And PLLA MS was easy to cause redness and swelling due to the strong inflammation after injection, the strong inflammation was not conducive to recover which was a major disadvantage of PLLA as a filler.

Collagen fibers are the most widely distributed fiber in connective tissue. In order to evaluate the effect of chiral PLA microspheres on collagen regeneration, collagen in the sections were observed under the optical microscope. Masson staining is an authoritative and classic technique for the staining of collagen fiber; it is frequently used in the analysis of collagen and smooth muscle. Collagen fiber, mucus, and cartilage were stained blue, and cytoplasm, muscle, cellulose, and glia were stained red. Nucleus was stained blue purple [25]. By using picrosirius red staining, different types of collagen fibers can be distinguished according to their specific forms and distinct colors under polarized light. For example, type I collagen fibers are densely arranged showing red and yellow colors, while type III collagen fibers exhibit threadlike structure in green color [26]. As shown in Figs. 58 and 59

(Supporting information), on three days and seven days after injection of microsphere, the collagen fibers were arranged neatly and densely, with no significant difference among the three groups. However, on Day 30, the collagen fibers in the control group were messy and sparse, and the collagen loss was very serious. On the contrary, the collagen fibers in the PLLA MS group still remained normal arrangement. Meanwhile, the collagen density of PDLLA MS and PDLLA MS group was significantly sparser than that of PLLA MS group. Nevertheless, the results of PDLLA MS group stimulation of collagen regeneration were better than that in PLLA MS group (Fig. 4A). Similar results were shown in optical photo of tissue sections stained with picrosirius red. Three days after injection, a large amount of type I collagen fibers (red and yellow) were observed in all the groups (Fig. 510 in Supporting information). On Day 7, collagen fibers turned darker and sparser. Specifically, type III collagen fiber was generated in PLLA MS group which was showing green in Fig. 511 (Supporting information). This difference was enlarged on Day 30, when new collagens were

(Fig. 4B). The content of type I collagen in the PLLA group was still 84.7%, and type III collagen was 15.1%. Total of type I and type III collagen were 1.4 times than the collagen in the PDLLA group, and 1.1 times than the collagen in the PDLLA group (Fig. 4C). Which

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